

KING EXAMINATION COUCH

- 27615 KING EXAMINATION COUCH - white
- 27616 KING EXAMINATION COUCH - blue
- 27617 KING EXAMINATION COUCH - black
- 27619 KING EXAMINATION COUCH - white with hole
- 27618 KING EXAMINATION COUCH - green with hole

Strong and stylish structure in oval tube 50x25 mm painted with metallic gray epossidic powder. Top in washable, fire proof, leatherette material class 1 IM with high density padding. Head section adjustable by mechanical system. Maximum load: 130 kg. Size: 190x60xh 80 cm

Delivered in kit form. Made in Italy.

- 27511 COUCH ROLL HOLDER
- 27978 FOOT STOOL - 1 step
- 27979 FOOT STOOL - 2 steps



27619



27978

27979



WIDTH: 60 cm

MAX LOAD 130 kg

27616

EXAMINATION COUCHES AND FOOT STOOLS

- 27401 EXAMINATION COUCH - chromed - black
- 27402 EXAMINATION COUCH - white painted - grey
- 27444 EXAMINATION COUCH - white painted - green

STOOLS: see matching stools page 75

Suitable for medical office or small surgeries. Easy to assemble, following enclosed instructions.

Comes complete with mattress.

Upholstery in two sections: one 130x61xh 5.5 cm fixed and a 50x61xh 5.5 cm backrest adjustable through different positions.

Delivered in kit form. Made in Italy.

- 27418 COUCH ROLL HOLDER

Adjustable roll holder suitable for round (Ø27-30 mm) or square (27-30 mm x 27-30 mm) metal couch leg.

- 27404 CLASSIC EXAMINATION COUCH - chromed - beige

- 27405 CLASSIC EXAMINATION COUCH - chromed - blue

- backrest with hole

- 27445 CLASSIC EXAMINATION COUCH - chromed - apricot

Chromed mainframe and legs constructed in chrome-plated steel. Legs fitted with plastic feet.

The two-section top has resilient foam padding covered with washable and fire-proof plastic material.

Backrest fully adjustable.

Available in cream or blue colour.

Other colours on request, minimum 10 pieces.

Size: 180x57xh 73 cm.


Maximum load 100 kg.

Delivered in kit form.

Made in Italy.

- 27408 MOUTH-NOSE HOLE PLUG for code 27405 - blue

STOOLS: see matching stools page 115



Couch size
180 x 61 x h 75 cm


27401

27402

Backrest adjustable through different position



Simulated leather, black/grey colour washable upholstery in two sections



27405

27404



MASSAGE BED



27421

- 27421 MASSAGE BED

Comfortable massage bed, ideal for physical therapy and wellness centres.

Back rest adjustable.

Supplied with couch roll holder.

Product size: 183x63x64 cm

Packing size: 188x17x66 cm

Net/Gross weight: 28/31 kg