

iHEALTH WIRELESS PROFESSIONAL HEALTH MANAGEMENT


9 SINGLE LEAD WEARABLE AMBULATORY ECG

• **23535 SINGLE LEAD ECG HOLTER** - available June 2014

Advised for use on patients who may be asymptomatic or who may suffer from transient symptoms such as palpitations, anxiety, syncope, light-headache, dizziness, fatigue, pre-syncope and shortness of breath.

It is simple to use, the patient launches the app on iPhone and it begins to record.

- designed to take continuous ECG measurements from the chest for 24 hour monitoring
- connects to a mobile device via Bluetooth 4.0 BLE to view the real-time ECG waveform
- saves data to your PC or mobile device
- internal memory stores up to 3 days of data
- thin and small ergonomic design to ensure comfort and good electrical contact, sticks directly to the skin
- rechargeable battery and sweat proof
- free iHealth mobile app and cloud service includes personalized health tools.



23535

CE 0197

TECHNICAL SPECIFICATIONS

Single Lead;
Power: 3.7 V Li-Ion, 200 mAh
Sample rate: 250 Hz
Bandwidth: 0.67 to 40 Hz
Input impedance: > 10 Mohm
CMRR > 60 dB
Heart rate range: 30 to 280 bpm
Size: 115x40x10 mm
Bluetooth 4.0 low energy

MOBILE DEVICE COMPATIBILITY

iOS 5.0
- iPhone 4S+
- iPad mini
- iPad 4+
- iPod touch (5th generation)

10 WEARABLE AMBULATORY BLOOD PRESSURE MONITOR

• **23540 WEARABLE AMBULATORY B.P. MONITOR** - available Sept. 2014

- non invasive method of obtaining blood pressure readings over 24 hours
- measures and records systolic/diastolic numbers and heart rate
- programmable blood pressure measurement cycle
- connects to a mobile device via dual bluetooth 3.0 & 4.0 BLE
- free iHealth mobile app and cloud service

Supplied with disposable arm cuff 22-42 cm

• **23541 DISPOSABLE ARM CUFF 22-42 cm** - spare

TECHNICAL SPECIFICATIONS

Rechargeable battery
Programmable measurement cycle setting (15, 30, 45 min, 1, 2 hours)

MOBILE DEVICE COMPATIBILITY

iOS 5.0
- iPhone 4+
- iPad mini
- iPad 2+
- iPod touch (4th generation)



23540

iBABY MONITOR FOR BABY WELLNESS - APPLE COMPATIBLE

iBABY APP - free



• **23480 WIRELESS BABY CAMERA iBABY M2**

The iBaby M2 Monitor follows your little one with a swipe of your iOS touchscreen.

See your baby through your mobile device and when baby's movement or cry activates alerts you can respond immediately.

The two-way audio feature lets you hear and speak to your little one with no external speaker wherever you are.

In addition to its unique and modern design, this monitor includes a built-in rechargeable battery which allows you to move the camera from room to room.

Bluetooth connectivity for iPad, iPod, iPhone.

Package and manual in 5 languages (GB, FR, IT, ES, DE).

2 years warranty.

BG1 BUTTON SIZE GLUCOSE MONITOR - MOBILE PHONE TRANSMISSION

• **23509 BG1 BUTTON SIZE GLUCOSE MONITOR**

Button size, economic glucose monitor with advanced technology and good reliability. A convenient way to measure your blood glucose levels and transmit by your phone.

Ear phone universal plug compatible with any mobile.

Free Apple and Android App: iHealth Glucosmart and iHealth MyVitals.

CE 0197



GIMA code	GLUCOSE MONITOR ACCESSORIES
23511	Glucose strips - box of 25
23512	Glucose strips - box of 50
23916	Lancets 28G - box of 100
23915	Auto lancets device "OnCall"

TECHNICAL SPECIFICATIONS

Size: 53x30x96 mm (2.1"x1.2"x0.38")
Measuring method: Amperometric technology using glucose oxidase
Result range: 20 mg/dL~600 mg/dL (1.1 mmol/L~33.3 mmol/L)
Power source: DC 3.0 V, CR1620
Blood source: Fresh capillary whole blood
Blood volume: Min. 0.7 micro liter